



SCYON Orthopaedics AG

International Implant card - Patient information

The implant card contains the information about the implants from SCYON used in a patient for the Total Hip Replacement (THR) procedure. The implant card is a credit card sized carton card that has a 2nd half folded behind the front. In total there are four sections.



Unfold the back portion to view the full card.




Front side, 2 sections
1 – General information
2 – Femoral stem information



Back side, 2 sections
3 – Femoral head information
4 – Acetabular cup information

Section 1

General information about the patient and implantation procedure. Positions 1, 2, and 3 will be completed by the healthcare provider. The 4th item is the web address where information about the manufacturer of the implants. Select the menu item “patient” to be directed to information targeted for patients of the THR procedure.



The image shows a form titled "SCYON Orthopaedics International Implant Card". It is divided into four numbered sections for completion:


- 1**: Patient identity, indicated by a person icon and a question mark.
- 2**: Date of implantation, indicated by a calendar icon with the number 31.
- 3**: Name and address of the implanting healthcare provider, indicated by a person icon with a plus sign.
- 4**: Information website, indicated by a laptop icon and the URL www.scyon.ch.

Item on card	Information	How completed
1	Identity of the patient	Handwritten by healthcare provider
2	Date of the implantation	Handwritten by healthcare provider
3	Name and address of the implanting healthcare provider	Handwritten by healthcare provider
4	Information website	Pre-printed on card


In Sections 2, 3, and 4 of the implant card, the patient labels provided with the implants are affixed to the card in the locations marked. The patient information card includes space for each of the three components required for the THR procedure from Scyon Orthopaedics—the femoral stem, femoral head, and acetabular cup. Therefore, one implant card covers all implants used.

Section 2

Information about the femoral component.



Section 2, Blank




Section 2, patient label affixed

Item on card	information	How completed
5	Device type in multiple languages	Pre-printed on card
6	Medical Device name	Pre-printed on label delivered with the implant
7	Unique Device Identifier (UDI) in human readable form	Pre-printed on label delivered with the implant
8	UDI code in AIDC format	Pre-printed on label delivered with the implant
9	Lot number	Pre-printed on label delivered with the implant
10	Name and address of the manufacturer of the implanted medical device.	Pre-printed on card

Section 3

Information about the femoral head component.



Section 3, Blank




Section 3, patient label affixed


Item on card	information	How completed
5	Device type in multiple languages	Pre-printed on card
6	Medical Device name	Pre-printed on label delivered with the implant
7	Unique Device Identifier (UDI) in human readable form	Pre-printed on label delivered with the implant
8	UDI code in AIDC format	Pre-printed on label delivered with the implant
9	Lot number	Pre-printed on label delivered with the implant

Section 4

Information about the acetabular component.



Section 4, Blank



Section 4, patient label affixed

Item on card	information	How completed
5	Device type in multiple languages	Pre-printed on card
6	Medical Device name	Pre-printed on label delivered with the implant
7	Unique Device Identifier (UDI) in human readable form	Pre-printed on label delivered with the implant
8	UDI code in AIDC format	Pre-printed on label delivered with the implant
9	Lot number	Pre-printed on label delivered with the implant
10	Name and address of the manufacturer of the implanted medical device.	Pre-printed on card

Appendix 1 – Symbol definitions

Symbol	Meaning
	Patient name or patient ID
	Information website for patients
	Unique device identifier
	Name and address of the implanting healthcare
	Date of implantation

Implant Card and Patient Information Symbols